



470267

IgE BioAssay™ ELISA Kit (Rabbit)

Kit Type

Sandwich ELISA

Tests 100

Sensitivity <0.8ng/ml

Detection Method

Colorimetric

Detection Range

3-400ng/ml

Sample Matrix

Serum, plasma, urine, cell culture

samples, biological fluids.

EU Commodity Code

38220000

Shipping Temp

Blue Ice

Storage Temp

4°C

Immunoglobulin E (IgE) is a type of antibody that has only been found in mammals. IgE is synthesised by plasma cells. IgE's main function is immunity to parasites. It has an essential role in type I hypersensitivity, such as allergic asthma, most types of sinusitis, allergic rhinitis, food allergies, and specific types of chronic urticaria and atopic dermatitis. IgE also plays a pivotal role in responses to allergens, such as: anaphylactic drugs, bee stings, and antigen preparations used in desensitization immunotherapy.

Sample Type

Serum, plasma, urine, cell culture samples, biological fluids.

Intended Use

This ELISA kit is used for in vitro quantitative determination of IgE.

Sensitivity

<0.8ng/ml

Range

3-400ng/ml

Specificity

Rabbit, no significant cross-reactivity or interference between this analyte and its

analogues was observed.

Precision

Intra-Assay: <10% Inter-Assay: <12%

Test Principle

This IgE BioAssay™ ELISA kit is a sandwich ELISA assay for the quantitative measurement of IgE in rabbit serum, plasma and cell culture supernatants in 90 minutes. The density of color is proportional to the amount of IgE captured from the samples.

Kit Components

Micro ELISA strip-plate Standard (360ng/ml), 1x500ul Standard diluent, 1x1.5ml HRP- Conjugate reagent, 1x6ml Sample diluent, 1x6ml Chromogen Solution A, 1x6ml Chromogen Solution B, 1x6ml Stop Solution, 1x6ml Wash buffer (30X),

1x20ml Plate sealers

Storage and Stability

Store powder at 4°C liquid at -20°C. Store other components at 4°C. Stable for 6 months For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap.

Intended for research use only. Not for use in human, therpeutic, or diagnostic applications.

References

No references available